

# OXYGEN CONCENTRATOR (Make in India Product)

## MO<sub>2</sub> TECHNICAL SPECIFICATIONS



### 10 LITRES

Model	MO2 (OCM 10A)
Voltage	230V~+/- 10%/ 50 Hz +/- 1 HZ
Highest flow	10 liter per min
Power consumption	770 VA
Weight	20.5 kgs
Oxygen output pressure	60 kPA - 80 kPA
Purity	90%~93%
Warranty	1+ 2 years extendable
Range of adjusting flow	0 to 10 liter per minute
Noise level	Less than 55 dB(A)
Dimension	(L* W* H) L 400* W 300* H 540 mm
Alarm	When 72% to 82%, the yellow indicator will shine, . When purity ≤ 72%, red indicator will alarm, display will show "Lo-O2". It will stop within one minute.

### 5 LITRES

Model	MO2 (OCM 05A)
Delivery rate	1~5 Litre Per Minute
Maximum Flow	5 litre per minute
Input Voltage	120 VAC (±10%)
Input Frequency	60 HZ
Highest flow	5 litre per min
Average power consumption	360 W
Power consumption	600 V A
Output pressure	58.6 ± 6 Kilopascals (kPa)
Weight	16.5 Kg
Purity	94%~95% at 2 Liter Per Minute 93%~95% at 3 Liter Per Minute 90%~93% at 5 Liter Per Minute
Noise level	< 46 dBA
Warranty	1+2 years Extendable
Dimension	L 400*W 300*H 510 (mm)
Alarm	High System Pressure, Low System Pressure, Power Failure, Low Oxygen Level, No Oxygen Flow,
Warranty	1+ 2 years extendable