



Provide Warmth
Anywhere, Everywhere

Embrace Nest is an easy-to-use, portable infant warmer that does not require continuous power and is used as supplement to KMC.

NICU



Transport



Maternity Ward



PHOENIX
Phoenix Medical Systems (P) Ltd.




Component Features

Embrace Nest Unique Benefits




WARMPAK™

Maintains temperature of ~37°C for at least 4 hours*
*duration of at least 4 hours at ambient temperature greater than 25°C
380 mm x 220 mm • 1.3kg



Special Phase Change Material in WarmPak™ maintains a temperature of ~37° C for at least 4 hours*



BABYWRAP™

Reusable and easy to sanitize using disinfectants or soap wash
520 mm x 250 mm • 0.3kg



Portable and simple to use




Only 15 Watts of power consumption in ready mode



ACCUTEMP HEATER

Precision heating mechanism prevents overheating of WarmPak™
440mm x 290 mm x 70 mm • 2.5kg • 230 V, 50 Hz




Does not require a constant supply of electricity



Easy to sanitize and reuse

Simple To Use



Product Order Information

Model no. (FD No.)	Product Name
FD0904	Embrace Nest Starter Pack
FD0903	Warmpak™ Refill Pack
FD0902	Accutemp Heater Refill Pack
FD0901	BabyWrap™ Refill Pack

Phoenix Medical Systems Private Limited

DP42, SIDCO Industrial Estate, Thirumudivakkam, Tamil Nadu – 600132, India
Toll Free India : 1800 425 3503, Phone: +91 44 2478 1705 / 2227,
✉ India : sales@pmsind.com, International: exports@pmsind.com
🌐 www.phoenixmedicalsystems.com