




BIOMIMICRY TECHNOLOGIES
PRIVATE LIMITED

**PULSEOM**TM



6302354783
www.biomimicrytechnologies.com



BIOMIMICRY TECHNOLOGIES
PRIVATE LIMITED



सिद्धिर्भवति कर्मजा





Specifications

Oxygen Level	: 35 ~ 99%
Pulse rate	: 30 ~ 220bpm
Perfusion Index	: 0.8 ~ 20%
Type	: OLED/LCD TFT Color Display
Parameters	: SpO ₂ , PR, PI, Waveform(Pleth,Line)

Working Condition

Temperature	: 5°C ~ 45°C
Humidity	: 15% ~ 80% non-condensing

Storage Condition

Temperature	: -10°C ~ 50°C (14°F - 122°F)
Humidity	: 10% ~ 90% non-condensing

6302354783
www.biomimicrytechnologies.com

