



SCHILLER

SOLUS[®] 1S

12.1" Multi-utility Patient Monitor



| NIBP | Manu | NIBP List | PR | Time |
|----------|------|---------------|----|-------|
| 160 / 90 | | 120 / 80 (93) | 60 | 14:20 |
| | | 120 / 80 (93) | 60 | 14:15 |
| | | 120 / 80 (93) | 60 | 14:10 |
| | | 120 / 80 (93) | 60 | 14:05 |

- ← Left
- Pat. Set
- NIBP
- Event
- Print
- Switch
- Freeze
- Review
- Right



FEATURES

- ❖ 12.1" high resolution colour TFT LCD display
- ❖ Total 8 waveform display
- ❖ ST and arrhythmia analysis
- ❖ 120 Hours of tabular and graphic trends & data storage
- ❖ Hemodynamic, Drug dosage, Oxygenation, Renal function & Ventilation calculation
- ❖ OxyCRG software
- ❖ Pacemaker detection facility
- ❖ SpO2 pulse-tone modulation (Pitch Tone)
- ❖ Rechargeable lithium-Ion battery
- ❖ Various mounting solutions
- ❖ Support HL7 protocols to HIS
- ❖ Central station connectivity : Wired



OPTIONAL

- ❖ Dual IBP
- ❖ EtCO₂ (Mainstream/ Sidestream)
- ❖ Thermal recorder
- ❖ Touchscreen
- ❖ CMS

Technical data are subject to change without notice

SCHILLER

The Art of Diagnostics

Manufacturing Unit - Unit III, AMTZ Campus, Ground Floor,
I HUB (B1) Building, Nadupuru Village, Pedagantyada Mandal,
Visakhapatnam, Andhra Pradesh - 530031
GSTIN : 37AADCS5091B2ZG

Manufactured at

AMTZ
सिद्धिभवति कर्मजा
www.amtz.in

