

GeneOx

TECHNICAL SPECIFICATIONS

Model	GeneOx
Weight	1.8 kg (approx)
Power Consumption	60 Watt (40W when not in use)
Rate Range	0 – 45 BPM (Titration mode)
Pressure Range	0 – 30 cm H2O
PEEP	Up to 30 cm
Battery backup	Up to 1 hour
Waveform	Flow / Pressure
Max leak flow compensation	10 L/min
Max flow rate	240 L/min
Mode- Non-invasive Mode	S / T / ST / BIPAP / CPAP
Supply Volatile	240 VAC
Supply Frequency	7 Hz – 55 Hz
Maximum Input Current	5.0 AMP
Internal Battery	11.1 V DC / 4.4 AH providing one hour backup

