

GeneOx

TECHNICAL SPECIFICATIONS

| Model | GeneOx |
|----------------------------|----------------------------------------------|
| Weight | 1.8 kg (approx) |
| Power Consumption | 60 Watt (40W when not in use) |
| Rate Range | 0 – 45 BPM (Titration mode) |
| Pressure Range | 0 – 30 cm H ₂ O |
| PEEP | Up to 30 cm |
| Battery backup | Up to 1 hour |
| Waveform | Flow / Pressure |
| Max leak flow compensation | 10 L/min |
| Max flow rate | 240 L/min |
| Mode- Non-invasive Mode | S / T / ST / BIPAP / CPAP |
| Supply Volatile | 240 VAC |
| Supply Frequency | 7 Hz – 55 Hz |
| Maximum Input Current | 5.0 AMP |
| Internal Battery | 11.1 V DC / 4.4 AH providing one hour backup |

