



**BIOMIMICRY TECHNOLOGIES**  
PRIVATE LIMITED

# Non-Contact Infrared Thermometer



Forehead  
Temperature




Rapid  
Measurement



Auto  
Shut Off



32 Memory



Object  
measurement



Fever Alarm

# Specifications

Measurement Object	Forehead, Surface, Room
Measurement Distance	1-10cm
Measurement Range	32-45C(89-120.2F)
Measurement accuracy	0.2-0.4(medical)
Resolution	0.1C

## Operating Environment

Relative Humidity	15-85%
Operating Temperature	5-45C

## Storage Environment

Storage Temperature	-20~ +55c
Relative Humidity	15~93%

